

Proton Therapy

*Se Byeong Lee, Ph.D.
Senior Medical Physicist*





NATIONAL CANCER CENTER
PROTON THERAPY CENTER




BASIC SPECIFICATIONS OF PROTON BEAM

History of Proton Therapy



- 1929 : Invented cyclotron, first particle accelerator, by Ernest Lawrence.
- 1946 : Proposed proton therapy by Robert R. Wilson.
- 1954 : The first proton patient (pituitary region, metastatic breast cancer??) at Berkley Laboratory.
- 1961, 1964, 1974 : Harvard Univ., U.C. Davis, Los Alamo tried proton therapy using cyclotron
- 1988 : Proton therapy got FDA approval.
- 1991 : The first hospital-based proton therapy facility using synchrotron in Loma Linda Univ.



Particle Therapy Center in the world



Electromagnetic energy loss of protons (Distal distribution)



2.2.2 The Bethe-Bloch Formula

The correct quantum-mechanical calculation was first performed by Bethe, Bloch and other authors. In the calculation the energy transfer is parametrized in terms of momentum transfer rather than the impact parameter. This, of course, is more realistic since the momentum transfer is a measurable quantity whereas the impact parameter is not. The formula obtained is then

$$-\frac{dE}{dx} = 2\pi N_a r_e^2 m_e c^2 \rho \frac{Z}{A} \frac{z^2}{\beta^2} \left[\ln \left(\frac{2m_e \gamma^2 v^2 W_{\max}}{I^2} \right) - 2\beta^2 \right]. \quad (2.26)$$

Equation (2.26) is commonly known as the *Bethe-Bloch formula* and is the basic expression used for energy loss calculations. In practice, however, two corrections are normally added: the *density effect* correction δ , and the *shell* correction C , so that

$$-\frac{dE}{dx} = 2\pi N_a r_e^2 m_e c^2 \rho \frac{Z}{A} \frac{z^2}{\beta^2} \left[\ln \left(\frac{2m_e \gamma^2 v^2 W_{\max}}{I^2} \right) - 2\beta^2 - \delta - 2 \frac{C}{Z} \right], \quad (2.27)$$

with

$$2\pi N_a r_e^2 m_e c^2 = 0.1535 \text{ MeVcm}^2/\text{g}$$

r_e :	classical electron radius = $2.817 \times 10^{-13} \text{ cm}$	ρ :	density of absorbing material
m_e :	electron mass	z :	charge of incident particle in units of e
N_a :	Avogadro's number = $6.022 \times 10^{23} \text{ mol}^{-1}$	β =	v/c of the incident particle
I :	mean excitation potential	γ =	$1/\sqrt{1 - \beta^2}$
Z :	atomic number of absorbing material	δ :	density correction
A :	atomic weight of absorbing material	C :	shell correction
		W_{\max} :	maximum energy transfer in a single collision.


The maximum energy transfer is that produced by a head-on or *knock-on* collision. For an incident particle of mass M , kinematics gives

$$W_{\max} = \frac{2m_e c^2 \eta^2}{1 + 2s \sqrt{1 + \eta^2 + s^2}}, \quad (2.28)$$

where $s = m_e/M$ and $\eta = \beta\gamma$. Moreover, if $M \gg m_e$, then


$$W_{\max} \approx 2m_e c^2 \eta^2.$$

Multiple Coulomb Scattering (Lateral distribution)





Calc. U. Weber, RKA AG

Nuclear interaction




Primary P loss ~ 1% per cm in water





Absorbed dose contributions in Proton therapy

- ✓ Incident Energy : 230 MeV
- ✓ Each line presents relative absorbed dose as a function of water depth.
- ✓ Dose contribution by alpha, d, t, is less than 1% of total dose.
- ✓ In the case of secondary particle, dose contribution of each slice is from 2% up to 12% before the bragg peak region.
- ✓ Absorbed dose beyond the bragg peak is mostly caused by secondary proton.




Relative Biological Effectiveness

(Recommend RBE(P) = 1.1 @ mid. Of SOBP)



ICRU 78




Dosimetry Tasks in Proton Therapy

- Acceptance testing and commissioning of treatment beam lines: *Vendor & User based on RFP.*
- Reference beam calibration of clinical beam
- Commissioning of treatment planning system: *Beam data measurements and TPS validation tests.*
- Periodic checks: *Daily, Monthly & Annual QA*
- Verification of dose delivery: *Patient QA*

PROTON BEAM DELIVERY NOZZLE SYSTEM


Schematic Plot of Beam Nozzle




- 1st scatter, 2nd scatter for **beam spread in lateral distribution**
- range modulator for **Spread Out Bragg Peak (SOBP)**
- ionization chamber 1,2 and Range verifier for **beam monitoring**
- snouts, aperture for **lateral beam shaping**
- compensator for **longitudinal beam shaping**



Beam Nozzle System (IBA universal nozzle)




Proton Beam modulation methods





Specifications of Modulated Proton Beam

Item	Singles Scattering	Double Scattering	Uniform Scanning	Pencil Beam Scanning
Range in Patient (g/cm ²)	3.35–20.4	4.51–28.42	3.42–32.1	0.5–32 (7.5–32)
Range mod. Adjustment, max. modulation (g/cm ²)	0.4(R>6), 0.17 9.2	0.2 20.75	0.5 Full	0.4 24
Range adjustment(g/cm ²)	0.09(R>6),0.05	0.1	0.1	0.4
Ave. dose rate (Gy/min)	5.93	3	1.15	~1.0/1L
Max. Field Size (cm)	4 (D)	24.3 (D)	40X30	22.5X25 At high E
Dose Uniformity (%)	1.25	1.05	1.5(R), 2.6(L)	3
Effective SAD (m)	2.55	2.19	2.12	2.12
Distal Penumbra (g/cm ²)	0.21	0.23	0.13	Beam sigma
Lateral Penumbra (cm)	0.17	0.48	0.30(y), 0.22(x)	4.9mm@R32cm 10.4mm@R8cm
	H&N, RS	General	Large Size	IMPT

RM + FS + IC set A




Second Scatterer System



Range Modulation wheels

Principle of SOBP




Beam Current Modulation (BCM) → Total 52 BCM records

Stop digit determines width of SOBP

Weighting_table stands for cyclotron current to determine flatness of SOBP

Range Modulator



#user commands

#Track (Large wheel rotation)

/GTR2/RM/track 5

#SOBP (Small wheel rotation)

1. position mode (pristine beam)

/GTR2/RM/angle deg 80.26


2. speed mode (SOBP)

/beam/bcm weighting_table SD


– SW 1 (Track 1,2,3)

– SW 2 (Track 4,5,6)

– SW 3 (Track 7,8,9)



SOBP on Scattering mode




4. Maximum Energy in Patient

2. Beam gating
signal
(10Hz)

5. RM wheel
synchronizing
(600 rpm)

1. BCM file pickup








Snouts
250, 180,
100 & Eye



Demonstration of Proton Beam Irradiation




PBS Beam Irradiation




IBA demonstration

Example of Patient Treatment



Tools for PBS measurement



3D Water Phantom



Lynx



Zebra(MLIC)






MatriXX



EBT Film


Concerning of Beam delivery Uncertainty



Ex.) 5mm beam off on Target (1 portal)
- Proton : 100% → 0% dose
- Photon : 100% → > 95% dose




Dose Cal. Error on TPS based X-ray CT



- Overall expected maximum error in SP_{rel} :
 - $\pm 1.8\%$ for bone
 - $\pm 1.1\%$ for soft tissue

	Soft tissue			Bone			Total
	Amount [cm]	wer [cm]	Abs. error [cm]	Amount [cm]	wer [cm]	Abs. error [cm]	Abs. error [cm]
Brain	10	10.3	0.11	1	1.8	0.03	0.14
Prostate	15	15.5	0.17	5	9	0.16	0.33

Patient Treatment Setup Verification



Example

Image
registration &
Position
Correction



- VARIAN A277 X-ray tube
 - kV range: 40 kV to 150 kV
 - Target angle: 7°
 - Focal spot sizes: 0.6 – 1 mm
 - mA/mAs range: 1 to 500 mAs
- VARIAN B150 housing
- VARIAN flat panels PAXSCAN 4030 R
 - Receptor type: Amorphous Silicon
 - Conversion screen: Gadolinium Oxyde (Lanex Fast)
 - Energy Range: 40kV- 150 kV
 - Active area: 29,3 cm x 40,6 cm
 - Pixels: 2304 x 3200
 - Resolution: 127 µm x 127 µm

Dosimetry Devices for PB



There are several tools, useful for scanning beam.

Some of them are essential and others can save your beam data measurement time.

Ex.) MLIC, Scintillation Plate type detector, 2D Array type Ionization Chamber, Large size PP type chamber, etc.

Maintenance PBT Facility

- **Biomedical Engineers : 9 man/year**
 - 1 for RF & PT system
 - 1 for Beam optics
 - 1.5 for Electronics
 - 1.5 for Mechanics
 - 1 for Software
 - 3 for Cyclotron OP
- **Medical Physicist**
 - 1 for Beam Calibration & Dosimetry support
- **IBA On-site Engineer**
 - 1 for General technical support

**** Preventive Maintenance : 320 items (including daily, weekly, monthly and yearly items)**


Weekly Schedule of PT Facility

time	Mon	Tue	Wed	Thu	Fri	Sat	Sun
6							
7							
8							
9							
10							
11							
12							
13							
14							
15							
16							
17							
18							
19							
20							
21	Patient QA & PM						
22							
23							
24							

Maintenance Statuses

- 9 Engineers & 3 shifts/day
- Machine : Proteus 235, Linac(trilogy, CL21EX, N21EX), SIMULATOR
- 10 years experience in PT


Recent Uptime per room(%)




Work Scope :

- 1. PT Facility maintenance**
 - CM(corrective maintenance) : ~75 times/y
 - PM(Preventive Maintenance)
 - : cyclotron PM 2 times(4days)/year,
 - Other part PM 2 times/year
 - Proton Accelerator operation
- 2. Linac-3EA, SIMULATOR maintenance**
CM : 50 times/year

Radiation Activity Distribution during Cyclotron PM






No.	Description	No.	Description
1	Opposite side of deflector	7	Rear end of gradient corrector
2	Front of Gradient corrector	8	Boundary of main coil back Gradient corrector
3	Right rear end of deflector	9	Boundary of main coil back permanent magnet
4	Boundary of main coil back deflector	10	Permanent magnet
5	Boundary of yoke back deflector	11	Boundary of yoke near permanent magnet
6	Left rear end of deflector	12	Front of fixed cavity tuner




- 6 times measurements(2014. 4~ 2017. 4)
- 1 day after beam shutdown(survey meter: 450B-DE-SI, VICTOREEN & FH40 G-10, THERMO)

Treatment Room Radiation Level



Radiation level monitor & area detector



Radiation Exposure statistics of radiation workers

Category	Exposure Dose in 2009
Therapists	< 1 mSv
Medical Physicists	< 1 mSv
MCR Operators	0.7 ~ 2.5 mSv
Maintenance Engineers	4.5 ~ 7.0 mSv



PROTON DOSIMETRY CODE OF PRACTICES


Words Definitions in Radiation Therapy

- **Iso-effective Dose:** Absorbed Dose*RBE, Cobalt-60 Equivalent Dose (CGy or GyE).
- **Absorbed Dose(Gy):** a physical quantity to measure the radiation energy absorbed by unit mass of substances. $1 \text{ Gy} = 1 \text{ J/kg}$.
- **Relative Biological Effectiveness(RBE):** the ratio of biological effectiveness of one type of ionizing radiation relative to Cobalt-60.

References



- ICRU Report 78, “Prescribing, Recording, and Reporting Proton-Beam Therapy”
- IAEA TRS no. 398, “Absorbed Dose Determination in External Beam Radiotherapy”.
- IAEA TRS no. 461, ”Relative Biological Effectiveness in Ion beam Therapy”.
- C. P. Karger, T. Kanai, et al, “Dosimetry for ion beam radiotherapy”, PMB 55(2010) R193-R234.
- C. Goma, et al, “Proton beam monitor chamber calibration”, PMB 59(2014) 4961-4971
- Ion Beam Application(IBA) Technical Reports.
- PTCOG Education Course Presentation.

International standards of D_w



Absorbed Dose Determination

- ICRU Proton Dosimetry code : Report 78
“Prescribing, Recording, and Reporting Proton-Beam Therapy”
- IAEA Proton Dosimetry code : TRS 398
“Absorbed Dose Determination in External Beam Radiotherapy: An international Code of Practice for Dosimetry based on Standards of Absorbed Dose to water”



Ionization Chamber Dosimetry protocol

- IAEA/ICRU recommending protocol based on Absorbed dose to water(D_w).
- Providing the formalism and the data to relate a calibration of an ionization chamber.
- Simple dosimetry system and easy to use in clinic.
- Availability to all kind of therapeutic beam modalities.

$N_{D,w}$ -based formalism

- Absorbed dose to water in a reference beam Quality : $Q_0(^{60}\text{Co})$
- Absorbed dose of Q(X-ray, Electron, Proton, Heavy ion).

$$D_{w,Q} = M_Q N_{D,w,Q_0} k_{Q,Q_0}$$

M_Q : Corrected Q Reading of Electrometer (PTP corr.).

N_{D,w,Q_0} : Calibration Coefficient of Ionization Chamber at ^{60}Co .

k_{Q,Q_0} : Beam Quality correction factor btw $Q_0 (^{60}\text{Co})$ & Q.

K_Q factor

- Measurements and MC Calculations based.
- I.C. perturbation of reference beam (Co) at PSSD/SSDL.

$$k_{Q,Q_0} = \frac{(s_{w,air})_Q}{(s_{w,air})_{Q_0}} \frac{(w_{air/e})_Q}{(W_{air/e})_{Q_0}} \frac{p_Q}{p_{Q_0}},$$

$$p = p_{cav} p_{cel} p_{dis} p_{wall},$$

P ~ 1 for proton
P = 1 for Carbon
P ≠ 1 for ⁶⁰Co

Table 4.11. Estimated relative standard uncertainties, u_c , of the calculated values for k_{Q,Q_0} for proton beams (IAEA, 2000).

Chamber type	Cylindrical		Plane parallel	
	Protons, u_c (%)	⁶⁰ Co + protons, u_c (%)	Protons, u_c (%)	⁶⁰ Co + protons, u_c (%)
$s_{w,air}$	1.0	1.1	1.0	1.1
Assignment of $s_{w,air}$ to beam quality	0.3	0.4	0.3	0.4
$W_{air/e}$	0.4	0.5	0.4	0.5
p_{cav}	0.3	0.3	0.3	0.3
p_{dis}	0.2	0.4	0.2	0.3
p_{wall}	0.6	0.8	0.6	1.6
p_{cel}	0.4	0.5	—	—
Combined relative standard uncertainty in k_{Q,Q_0}	—	1.7	—	2.1





$S_{w,\text{air}}$ & $W_{\text{air}/e}$ of Proton Beam

$S_{w,\text{air}}$ at ICRU 78



- ICRU 59 values
- PETRA MC calculations
(Spencer Attix cavity theory)
- Secondary electrons and
Nuclear inelastic processes.

*TRS 398 use the constant
value of 1.13 for proton.

$W_{\text{air}/e}$ at ICRU 78 & TRS 398 Inaccuracy reduction from Calorimeter and I.C. measurements



K_Q values for Ionization Chambers



$$R_{\text{res}} = R_p - Z_{\text{ref}}$$

Fig 10.2. Calculated values of k_Q for various cylindrical and plane-parallel ionization chambers commonly used for reference dosimetry, as a function of proton beam quality Q (R_{res}). Data from Table 10.III.

Reference Conditions


Table 4.7. Reference conditions for the determination of absorbed dose in proton beams (IAEA, 2000).

Quantity	Reference value or reference characteristics
Phantom material	Water
Chamber type	For $R_{\text{res}} \geq 0.5 \text{ g cm}^{-2}$, cylindrical and plane parallel chambers For $R_{\text{res}} < 0.5 \text{ g cm}^{-2}$, plane parallel chambers
Measurement depth, z_{ref}	Middle of the SOBP ^a
Reference point of the chamber	For plane-parallel chambers, on the inner surface of the window at its centres For cylindrical chambers, on the central axis at the center of the cavity volume
Position of the reference point of the chamber	For plane-parallel and cylindrical chambers, at the point of measurement depth z_{ref}
SSD	Clinical treatment distance
Field size at the phantom surface	$10 \times 10 \text{ cm}^2$, or that used for normalization of the output factors whichever is larger For small field applications (<i>i.e.</i> , eye treatments), $10 \times 10 \text{ cm}^2$ or the largest field clinically available



Recombination Corrections for proton

- Recombination (k_s) for **continuous** and **scanned continuous** beam
 - Two voltage method: **(V_1 / V_2)² instead of (V_1 / V_2)**
 - Energy / depth dependent



Reference dosimetry in a scanned pulsed proton beam using ionisation chambers and a Faraday cup


Stefan Lorin¹, Erik Grusell¹, Nina Tilly¹, Joakim Medin²,
Peter Kimstrand¹ and Bengt Glimelius¹

¹ Department of Oncology, Radiology and Clinical Immunology, Uppsala University,
Akademiska sjukhuset, SE-751 85 Uppsala, Sweden

² Department of Radiation Oncology, Copenhagen University Hospital, DK-2100 Copenhagen,
Denmark


*Independent Verification of
Recombination correction
needed by user!*

Reference Conditions for Spot Scanning



**Calibration at
SOBP**

Passive Scattering Beam



**Calibration at
SOBP**
*(or Calibration at
plateau)**


Spot Scanning Beam

* Reference condition are facility specific.
Required for TPS Beam data.

Proton spot scanning Calibration

PSI


(Pedroni et al, 2005)



MedAustron

$$D_w = (N/A) (S/\rho)_w * 1.602 \times 10^{-10}$$


Number of protons/MU




$$K(E) = \frac{N}{\text{MU}} = \frac{D_{\text{meas}}}{S_E(x)\text{MU}} \Delta x \Delta y$$

MD ACC

(Gillin et al, 2010)




Dose/MU



Measuring depth dose curve of Pencil beam

Detector size comparing to beam size matters!

- Bragg Peak Chamber is not large enough to collect signal from all protons present in the beam.



Correction parameter needed !

$$D_{\text{exp}}(R_0(k), z_p) = \text{corr}(R_0(k), z_p)^* D_{\text{mea}}(R_0(k), z_p)$$

$$\text{corr}(R_0(k), z_p) = \int_0^{\infty} r D_{\text{mea}}(r, z_p) dr / \int_0^{4.08} r D_{\text{mea}}(r, z_p) dr$$

Proton Beam Monitor Chamber Calibration




Figure 3. Ratio between the BMC calibration obtained with different ionization chambers—Farmer (□) and Advanced Markus (\diamond)—and the BMC calibration obtained with the Faraday cup. For the Farmer chamber, the BMC calibration is performed with both the reference point of the chamber positioned at z_{ref} (----) and the effective point of measurement of the chamber positioned at z_{ref} (—).

For the reference **Dosimetry of pristine bragg-peak proton beam, a plane-parallel I.C. is recommending rather than a cylindrical I.C.**

Comparison of ICRU 59 and TRS 398

Table 4.14. Comparison of ICRU 59 (ICRU, 1998) dosimetry protocol and TRS 398 (IAEA, 2000) dosimetry code of practice for clinical proton beams.

Feature/quantity	ICRU 59 (ICRU, 1998)	TRS 398 (IAEA, 2000)
Ionization chamber	Cylindrical	Cylindrical and plane-parallel ($R_{\text{res}} \geq 0.5 \text{ g cm}^{-2}$) Plane-parallel ($R_{\text{res}} < 0.5 \text{ g cm}^{-2}$)
Wall material	Graphite or A-150 plastic	Graphite for cylindrical chambers
Gas filling	Ambient air	Ambient air
Chamber volume	>0.5 cm ³ for beams $\geq 5 \text{ cm}$ diameter ~0.1 cm ³ for beams $< 5 \text{ cm}$ diameter	–
Water proof sleeve	–	$\leq 1 \text{ mm}$ thick PMMA
Dose specification	Water	Water
Calibration beam	^{60}Co	^{60}Co
Calibration coefficient	Primarily N_K , also N_X , $N_{D,w}$	$N_{D,w}$ only
Beam quality	Residual range (to distal 10 % level)	Residual range (to distal 10 % level)
Phantom material	Water (or other materials which match electron density of water)	Water
Reference point for measurement	Middle of SOBP	Middle of SOBP
Field size	–	Depth of 3 g cm^{-2} for plateau irradiations
SSD	–	$10 \times 10 \text{ cm}^2$
Stopping powers	ICRU (ICRU, 1993a)	Clinical treatment distance
$(w_{\text{air}}/e)_p$ (J C ⁻¹)	34.8 ± 0.7 (ambient air)	PETRA (Medin and Andreo, 1997b)
$(w_{\text{air}}/e)_{60_{\text{Co}}}$ (J C ⁻¹)	33.77 ± 0.05 (ambient air)	34.23 ± 0.13 (dry air)
$(s_{w,\text{air}})_{60_{\text{Co}}}$	1.134	33.97 ± 0.07 (dry air)
Chamber perturbation factors	No	1.133
Relative uncertainty in absorbed dose (1σ)	2.6 % (Jones, 2001d)	Yes 2.0 % (cylindrical chambers) 2.3 % (plane-parallel chambers)

Practical absolute dosimetry formula

- Measurement in water phantom.
- Detector of choice:
 - PPC ($R_{res} < 0.5 \text{ g/cm}^2$)
 - CC ($V_s \geq 0.6 \text{ cm}^3$) ($R_{res} > 0.5 \text{ g/cm}^2$)
 - Large \mathbf{V} : $\uparrow V \rightarrow \uparrow \text{Sensitivity} \rightarrow \uparrow \text{SNR}$



$$D_w = M \times N_{Dw} \times k_q \times k_{tp} \times k_s \times k_{pol} \times k_{el}$$

M – collected charge

N_{Dw} – calibration factor (*certificate*)

k_q – beam quality

k_{tp} – t/p correction factor

k_s – recombination factor

k_{pol} – polarity factor

k_{el} - calibration factor (*electrometer*)

$$k_{tp} = \frac{T \times p_0}{T_0 \times p}$$

$$k_s = a_0 + a_1 \left(\frac{M_1}{M_2} \right) + a_2 \left(\frac{M_1}{M_2} \right)^2$$

$$k_{pol} = \frac{|M_+| + |M_-|}{2M}$$

Standard condition

$$p_0 = 101.325 \text{ kPa}$$

$$t_0 = 20^\circ\text{C} (\text{in EU})$$

$$t_0 = 22^\circ\text{C} (\text{in US})$$

$$T = 273.15 + t$$

Worksheet for Dw in Proton beam

**Worksheet for the determination of the absorbed dose to water
in a proton beam**

User: [REDACTED] Date: [REDACTED]

1. Radiation treatment unit and reference conditions for $D_{w,Q}$ determination

Proton therapy unit: **Gantry2 (Famer chamber)**

Nominal dose rate: **200.0 MU min⁻¹** Nominal energy: **15.0 MeV**
 Reference phantom: **water** Practical range, R_P : **5 g cm⁻²**
 Reference field size: **10 cm x cm** Width of the SOBP: **5 g cm⁻²**
 Reference depth z_{ref} : **12.5 g cm⁻²** Reference SSD: **cm**
 Beam quality, $Q(R_{res})$: **2.5 g cm⁻²**

2. Ionization chamber and electrometer

Ion. chamber model: **PTW 30006 / 30013** Serial No.: **[REDACTED]** pp cyl
 Chamber wall material: **PMMA** thickness: **0.057 g cm⁻²**
 Waterproof sleeve material: **[REDACTED]** thickness: **[REDACTED] g cm⁻²**
 Phantom window material: **[REDACTED]** thickness: **[REDACTED] g cm⁻²**

Absorbed-dose-to-water calibration factor $N_{D,w}$: **0.05398** Gy/nC Gy/rdg
 Calibration quality Q_0 : **Co-60**

Reference conditions for calibration
 P_0 : **101.3 kPa** T_0 : **22.0 °C** Rel. humidity: **50 %**
 Polarizing potential V_1 : **400 V** Calibration polarity: +ve -ve corrected for polarity effect
 User polarity: +ve -ve

Calibration laboratory: **[REDACTED]** Date: **[REDACTED]**
 Electrometer model: **[REDACTED]** Serial no.: **[REDACTED]**
 Calib. separately from chamber: yes no Range setting: **[REDACTED]**
 If yes Calibration laboratory: **[REDACTED]** Date: **[REDACTED]**

3. Dosimetry reading ^a and correction for influence quantities

Uncorrected dosimeter reading at V_1 and user polarity: **45.3** nc rdg
 Corresponding accelerator monitor units: **200** MU
 Ratio of dosimeter reading and monitor units: $M_1 = 0.2265$ nc/MU rdg/MU

(i) P: **101.2 kPa** T: **22.5 °C** Rel. humidity: **50 %**
 $k_{T,P} = \frac{(273.2 + T)}{(273.2 + T_0)} \frac{P_0}{P} = 1.003$

(ii) Electrometer calibration factor k_{elec} : **[REDACTED]** nC/rdg dimensionless

(iii) Polarity correction ^b rdg at $+V_1$: $M_+ = 40.65$ rdg at $-V_1$: $M_- = 40.62$

$$k_{pol} = \frac{|M_+| + |M_-|}{2M} = 1.000$$

(iv) Recombination correction (two-voltage method)
 Polarizing voltages: V_1 (normal) = **400 V** V_2 (reduced) = **200 V**
 Readings at each \circ V: $M_1 = 40.65$ $M_2 = 40.58$
 Beam type: pulsed pulsed-scanned
 Voltage ratio $V_1 / V_2 = 2.0000$ Read. ratio $M_1 / M_2 = 1.002$
 $a_0 = 4.7110$ $a_1 = -8.2420$ $a_2 = 4.5330$

$$k_s = a_0 + a_1 \left(\frac{M_1}{M_2} \right) + a_2 \left(\frac{M_1}{M_2} \right)^2 = 1.003$$

Corrected dosimeter reading at the voltage V_1 :
 $M_Q = M_1 k_{TP} k_{elec} k_{pol} k_s = 2.2786E-01$ nC / MU rdg / MU

4. Absorbed dose rate to water at the reference depth, z_{ref}

Beam quality correction factor for user quality Q : $k_Q = 1.0310$
 taken from Table 31 Other, specify: **[REDACTED]**

$$D_{w,Q}(z_{ref}) = M_Q N_{D,w} k_Q = 1.2681E-02 \text{ Gy / MU}$$

Basic Beam Parameters

- Relative:
 - Depth dose curve.
 - Beam profile.
→ Film, Scintillating screen, Ionization chamber
- Absolute:
 - Beam output.
→ Ionization Chamber, Calorimeter


Clinical Parameters

- 100% dose level : average dose in uniformly treated region
- Range : 90% distal dose point
- Field size : width of 50%-50% dose, treatment field size
- Lateral penumbra : width of 80%-20% lateral dose
- Distal fall-off : width of 80%-20% distal dose
- Dose uniformity of Single-field Uniform Dose
- Beam spot asymmetry :


$$\left(\frac{(\sigma_x - \sigma_y)}{(\sigma_x + \sigma_y)} \right)$$

Field Definition

SOBP (90%–90%)



Field Size (50%–50%)






Summary

- **Code of Practice TRS 398 (IAEA 2000, N_{dw} –base formalism) is widely adopted for proton beam dosimetry in the world.**
- **Still, the dosimetry standard of spot scanning proton beam is not sufficient.**
- **For proton spot scanning and IMPT trial, there are lack in the specialized dosimetry tool.**
- **Proton dosimetry has a dependency of beam delivery system and home developed practices are needed.**



DOSIMETRY TOOLS FOR PROTON THERAPY


Detectors for Proton Beam Measurement




Detectors for absorbed dose applied in ion beam radiotherapy*

Detector	Advantage	Disadvantage	Application in ion RT
Calorimeter	Direct dose measurement	High effort, knowledge of chemical heat defect and thermal heat conduction required	Potential primary standard in the future, k_Q measurements
Ionization Chamber	High accuracy and reproducibility, small LET and energy dependence, easy to handle, many chamber types for different Applications	deviation from calibration conditions required, incomplete knowledge of corrections (chamber dependent)	Reference dosimetry, commissioning, dosimetric QA, dose verification, beam Monitoring
Radiochromic films	High spatial resolution, 2D measurement, Linear response, no daylight sensitivity, self-developing, less LET- and energy Dependent	LET and energy dependence, dose cannot be obtained from optical density in mixed fields, off-line analysis required Complex evaluation protocols, long term self-development, mechanical sensitivity	Measurement of lateral dose profiles, beam widths, field geometry and homogeneity, documentation of beam ports
Silicon diode and diamond	High spatial resolution, high signal, electronic read-out	LET, dose rate and energy dependence	Lateral profile Measurements
TLD	High spatial resolution	LET and energy dependence, off-line Evaluation	Point measurements, <i>in vivo</i> dosimetry
OSL detector	High spatial resolution, linear response, repeated electronic read-out	LET and energy dependence	Point measurements, Profiles
Scintillating screen	High 2D spatial resolution, linear intensity-independent response, electronic read-out	LET and energy dependence, large device	1D/2D distributions, field homogeneity, beam width, dose Verification

* "Dosimetry for ion beam radiotherapy", C.P. Karger, ... T. Kanai, P.M.B. 55(2010)R193-R234




Depth Dose Measurement at Small field beam



Careful selection of detect for small beam size needed.

Scanning Beam: depth dose curve

- Water phantom with a **single detector** is extremely slow in case of scanning beam modes: Full modulation has to be rerun for each ROI.
- Multi-Layer Ionization Chamber(**MLIC**): whole DD should be captured at once.



Beam Profile Measurement


- 2D Dose measurements for QA / validation purposes.

MatriXX



1020 ICs covering area of 24.4x24.4 scm
Distance between ICs = 7.6 mm
I.C. Size = 4.5 mm, I.C. Height = 5 mm

Not recommended for high-res measurements:
evaluation of penumbras
spot size measurements
coincidence measurements



Beam Profile Measurement


- 2D Dose measurements for QA / validation purposes.

Lynx

Scintillating screen + CCD camera

Good for high-res measurements: Evaluation of penumbras,
Spot size measurements, Coincidence measurements

But.. Exposure control needed , Large detector size,
Relative dose only, Quenching effect exist.





Beam Profile Measurement

- 2D Dose measurements for QA / validation purposes.


Self developing Film, Easy to use,
High resolution 2D dos distribution,
Quenching Effect, Relative dosimetry,
Field size measurement,

GafChromic
Film





PATIENT QA

Patient dedicate Device QA



Aperture / Range Compensator



Mechanical checks are applied

1. opening shape overlap on the plan sheet.
2. RC x-y coordinate & thickness measurement/comparison
3. Rotation and translation checks


Treatment Field QA

- **Passive Beam Plan:**
 - Range/SOBP(Optional): fine beam parameter finding, Markus chamber & water phantom, MLIC, Peak finder.
 - Output: Verification Plan, Farmer chamber & water(or solid water) phantom, without range compensator.
 - 2D dosimetry(if needed): Film, water phantom, MatriXX.
 - In vivo dosimetry: Film(EBT).
- **Spot Scanning Plan:**
 - MU confirmation: Verification Plan, Ion chamber/MatriXX & solid water phantom, 2 or more different depths.
 - 2D dosimetry: Verification Plan, MatriXX, solid water phantom, 2 or more different depths.
 - Depth dose: Central region with Markus/MatriXX and Solid water phantom.

2nd check method : Output factor, Range, SOBP calculation


Output

Hanne Kooy's Paper



Ψ : Output factor
 Ψ_c : Output factor (calib)
CF : Constant related
 RM opt.
D_{0,c} : Entrance Dose (calib)
r : $r = (R-SOBP) / SOBP$
a₀, a₁ : theoretical values


Range



Y : Machine nominal range
X : measured range
 1D linear Fitting func.

$$\rightarrow Y = P1 \cdot X + P2$$

SOBP



Y : Mach. SOBP – Measur. SOBP
X : Measur. SOBP
 4D polynomial fitting func.

$$\rightarrow Y = P1 \cdot X^4 + P2 \cdot X^3 + P3 \cdot X^2 + P4 \cdot X + P5$$

* It could be also used instead of beam measurement procedure.

Passive beam P-QA example:

Range/SOBP fine tuning finding & Output Measurement

1. RTP parameter

	Plan	Verification	Dose dist.
Prescription Dose :	235 CcGE		
Range in patient :	25.05 g/cm ²		25.07 g/cm ²
SOBP width :	10.39 g/cm ²		10.27 g/cm ²
Snout position :	27.9 cm	27.9 cm	
Air Gap :	2 cm	12 cm	
SCD of V.P. :	230.6 cm	234.4 cm	
Dose at V.P. :	241.4 CcGE	230.8 CcGE	

2 . Range/Mod calibration

	Convalgo		Measurement			
	Rp	SOBP	Rp	SOBP	dRp	dSOBP
1	25.05	10.39	25.10	10.27	0.09	0
2	24.96	10.39	25.07	10.3	0	0.03
3						
4						
5						
final	24.96	10.39	25.07	10.3	0	0.03



3. OutPut Calibration

Range in patient :	25.07 g/cm ²	RBE Correction :	1.1
SOBP width :	10.3 g/cm ²	Physical Dose :	209.8 cGy
Reference depth :	20.4 g/cm ²	Output factor :	1.4164 cGy/MU
Practical Range :	25.88 g/cm ²		
Calibration MU :	200 MU		
Fn :	1.045927	▷ Treatment MU :	148.13 MU

4. Summary - Field setup parameter

Range in patient :	24.96 g/cm ²	Field radius :	B8 7 cm
SOBP width :	10.39 g/cm ²	Gantry Angle :	270 °
Treatment MU :	148.1 MU	Snout extension :	27.9 cm
Dose rate :	2 Gy/min	Tolerance table :	TOL8

5. Depth Dose Curve of calibrated field



Measurement: line

NCC, Korea example

Spot beam Scanning P-QA example:

2D dose analysis, point dose comparison, etc.




matriXX & S.W.P

Large dose gradient and small volume distribution.

→ difficult to make a measurement with current dosimetry tools that have volumes and setup uncertainties.

Measured dose plane



Dose profile comparison

Calculated dose plane

γ - analysis result

of Fields : 2
Modulation Technique : MFO

1A LAO G30,CO RS74			
Actual Depth (cm)	2.0	5.0	8.0
Max dose (cGy)	130.67	128.02	-120.49
Init Gamma	91.59	96.74	100.00
AVG. Diff.	2.60	1.10	0.20
Final Renorm.		-2.00%	
Final Gamma	97.32	99.04	100.00
Final Avg. Diff.	1.18	0.13	-0.42
Renorm. 1% DOWN	95.41	97.93	100.00
Renorm. 2% DOWN	97.32	99.04	100.00
Renorm. 3% DOWN			

Temperature : 22.60 °C
Pressure : 1006.30 hPa

Final Renormalization	Plan Normalization	MU
Before correction (original plan)	99.100	429.17
After correction	101.122	420.59
Calculated Value (plan per field)	101.122	420.59
Final Plan Value (plan per field)	101.12	420.51

1B LPO G100,CO RS74			
Actual Depth (cm)	2.0	5.0	8.0
Max dose (cGy)	154.17	174.77	157.90
Init Gamma	97.32	97.66	99.74
AVG. Diff.	1.30	1.20	-1.50
Final Renorm.		-2.00%	
Final Gamma	99.99	100.00	97.00
Final Avg. Diff.	0.16	0.01	-1.61
Renorm. 1% DOWN	98.78	99.13	100.00
Renorm. 2% DOWN	99.99	100.00	97.00
Renorm. 3% DOWN			

Final Renormalization Plan Normalization MU

Final Renormalization	Plan Normalization	MU
Before correction (original plan)	99.100	568.53
After correction	101.122	557.16
Calculated Value (plan per field)	101.122	557.16
Final Plan Value (plan per field)	101.12	557.42



PROTON THERAPY FACILITY IN NCC

NCC Project schedule

- Release of final RFP : 2002. 04
- Contract with IBA : 2002. 06
- Ground breaking : 2003. 03. 21
- Equipments installation : 2004. 11 – 2005. 06
- Beam Production : 2005. 10. 21
- Beam tuning (Cyclotron → ESS → BTS → TR)
: 2005. 09 – 2006. 07
- Treatment Beam Calibration (GTR1→FBRT1→GTR2)
: 2006. 05 – 2007. 04

경

국민을
보호하는
로

축

국민의 알맞은
암환자의 삶의
로


국립암센터 연구동 및 양성자 치료시설 기공

2003. 3. 21 (금)



NCC Project History

- Release of final RFP : 2002. 04
- Contract with IBA : 2002. 06
- Ground breaking : 2003. 03. 21
- Equipments installation : 2004. 11 – 2005. 06
- Beam Production : 2005. 10. 21
- Beam tuning (Cyclotron → ESS → BTS → TR)
: 2005. 09 – 2006. 07
- Treatment Beam Calibration (GTR1→FBRT1→GTR2)
: 2006. 05 – 2007. 04



NCC Project History




- Acceptance test of Gantry treatment Room 2 : 2006. 10
- Acceptance test of Fixed beam treatment Room 1 : 2006. 12
- Acceptance test of Gantry treatment Room 3 : 2006. 12
- GTR1 (commissioning & Data taking) : 2006. 11 – 2007. 03
- 1st patient in GTR2 : 2007. 03. 19
- 1st patient in GTR3 : 2007. 09
- 1st patient in Fixed Beam Room : 2009. 06

NCC Project History

- Start of International Patient Treatment: 2010. 11
- National insurance reimburse the proton therapy
for pediatric cancer (< age 18): 2011. 04
- National Insurance reimburse the proton therapy
for most of cancers : 2015. 09
- 1st Pencil Beam Scanning Treatment in GTR3: 2015. 09
- 2400 patient treated in Proton Treatment Rooms by 2017

Statistics of PT Patients

(2007. 3 ~ 2017. 12 / 2,478 patients)



Proton therapy Facility



Proteus 235 Specification

Weight : 220 ton Height : 210 cm Diameter : 434 cm

Energy : 230MeV

Max. extracted beam current : 300nA

RF frequency : 106 MHz



Cyclotron



Proteus 235 Specification

Weight : 220 ton Height : 210 cm Diameter : 434 cm

Energy : 230MeV^{PM}

Max. extracted beam current : 300nA



RF frequency : 106 MHz

Energy Selection System





Transmission efficiency ESS



Max Beam current in treatment room = 300 nA x Transmission ESS


* transmission ESS to treatment room = ~100%

* $300 \text{ nA} = 4.8 \times 10^{12} \text{ protons/sec.}$


Beam Transportation System



Gantry Treatment Room - 1



Gantry Treatment Room - 2



Fixed Beam Treatment Room



Experiment Room



Main Control Room : operation consoles



MCR - Beam Production

Cooling
Vacuum
Magnet
RF
Radial Probe
Ion Source
ISEU
Deflector
Extraction Quads
Extraction Steering

PLC <-> P.S. Serial Link Fail

ACU Health status: NOK
ACU counter: 0

Cyclone 235 General Status

Main Coil PS	I Main Coil: 0.000 A	Source PS	Arc Tuning	I Fil.: 0.0 A
ON OFF	I Coil 1-3: 0.000 A	ON OFF	V Arc: 0.00 V	I Arc: 0.00 mA
	I Coil 2-4: 0.000 A			
Main Coil Tuning	I Compensation Coil: 0.00 A	Motor Position:	0 Digits	

LLRF → SSA → IMPA → FPA → RF Cavity

External Interlocks X External Conditions X External Conditions

ON Sby OFF

Forw Pw: 0.00 kW
V Dee 1: 0.00 kV
V Dee 2: 0.00 kV

Vdee Setpoint: 0 Digits

Deflector PS

V PS: 0.0 kV	Radial Probe	Position: 0.0 mm
I PS: 0.0 mA		Energy: 0.0 MeV
ON OFF	Motor1 Pos.: 0 Digits	In Limit Switch
	Motor2 Pos.: 0 Digits	Out Limit Switch

q1&q2

I Q1: 0.0 A
I Q2: 0.0 A

ON OFF

ALARM
Safety Not OK

Steering 1 PS

ON OFF

Steering 2 PS

ON OFF

Save Eqt Parameters

Edit Eqt Parameters

ESS layout

ESS Magnets Control

Cyclo Vacuum

Penning1: 0.0E+00 mbar

Hydrolic Station


Cooling

Instan

Main Control Room : ISEU, RF, Magnet E



Building Interlock



Monitoring & Safety System Interlock

Safety System Interlocks Summary

Computer Health

TSS	SCU	TGS	ACU
ECU			BLPSCU
SCU			
MCRS1			
TCU1			
PCU1			
TCU2			
PCU2			

ALL HEARTBEATS PRESENT

Room Security

TSS	SCU	SubZone 2	Cyc Vault (SZ1)
			MCR Master Key
		SubZone 3	
		SubZone 4	
		TCR1 Master Key	
		TR 1	
Beam Ready Time-out elapsed			
		Gantry 2	
		TCR2 Master Key	
		TR 2	
Beam Ready Time-out elapsed			
		Open/OFF	Secure/ON

Accelerator

R. Probe Center
Beam Current
Redund. Shutoff
RF Enabled
ESS BS2 Enabled
IS Enabled

Crash Buttons

MCR	SZ1	SZ2	TR1/ TCR1	SZ3	SZ4	TR2/ TCR2	GTR2
G	G	G	L	G	G	L	G
G = Global effect				L = Local effect			

Gantry

TR1	TR2	Crash Status
		Access Pt PreC
		Snout locked
		CW LS
		CCW LS
		Cab1 Press.
		Cab2 Press.
		Brk Press.
		PCU Sgle PPD Sel
		PCU Brk Rel Rq
		Drv1 Enabled
		Drv2 Enabled
		Motion Enable

PPS

TR1	TR2	PGU Sgle PPD Sel
		Access Pt PreC
		Crash Status
		Collision Memory
		Motion Enable

Snout

TR1	TR2	PCU Sgle PPD Sel
		Access Pt PreC
		Crash Status
		Motion Enable

Treatment

TR1	TR2	XREU
		ICEU1
		ICEU2
		ICEU3
		DCEU
		FSEU
		SSEU
		VCEU
		RVEU
		STEU/SREU
		RMEU (scatt.)
		SMEU (scatt.)
		SMPS OFF
		RMEU (wobbl.)
		SMEU (wobbl.)
		SMPS OK
		Service Mode ON
		TR NOT Selected

DIDs (Flat Panels)

TR1	TR2	PCU Sgle PPD Sel
		Access Pt PreC
		Crash Status
		Motion Enable

X-Rays


TR1	TR2	Ax 1 tube in pl.
		Ax 1 IP in pl.
		Ax 2 IP in pl.
		Crash Status
		Axis 1 Enable
		Axis 2 Enable

**SOME CIRCUIT BREAKER(S) TRIIPPED IN THE TSS.
THE INFORMATION SHOWN HERE MAY NOT BE ACCURATE !**

Detailed Troubleshooting

TR2 GANTRY HARDWARE & OTHER CONDITIONS		TSS	SCU
	✓ PCU health status OK Press the hs.reset button in TSS	✓	✓
	✓ SCU health status OK Press the hs.reset button in TSS	✓	✓
	✓ MCRS1 health status OK Press the hs.reset button in TSS	✓	✓
	✓ Snout locked	✗	✓
	✓ CW limit switch		✓
	✓ CCW limit switch		✓
	✓ Brakes pneumatic cabinet 1 supply pressure		✓
	✓ Brakes pneumatic cabinet 2 supply pressure		✓
	✓ Brakes air pressure		✓
NOT APPLICABLE (NO ACCESS POINT ACTIVE)		✓ Single patient positioning device selected in PCU	✓
		✓ Brakes release request from PCU	✓
		✓ Motor driver 1 enabled	✓
		✓ Motor driver 2 enabled	✓

Structure of the relay circuitry



Main Control Room : TSS relay racks




CCD camera system



CCD camera system



Radiation monitor with Area detectors



Milling Machine Shop



Thank you for your attention!

